

Gastric Lavage Tum-E-Vac®

Tum-E-Vac® Gastric Lavage Kits are fully assembled, disposable kits that promote quick, efficient lavage of stomach contents. This completely closed system is designed to protect healthcare professionals from potential exposure to stomach contents and charcoal solutions during gastric lavage or decontamination procedures. Available with or without charcoal. Latex free, single use, non-sterile.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Tum-E-Vac® Gastric Lavage Kit

- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Closed system helps prevent exposure to stomach contents and charcoal solutions
- Available with activated charcoal suspension (with or without Sorbitol)
- Twist-Lok™ expressing key aids in charcoal administration
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

Without Charcoal

Cat#	Description	Box
2054	Tum-E-Vac + 3L Solution Bag + 32Fr Lavage Tubing	10
2055	Tum-E-Vac + 3L Solution Bag	10
2075	Tum-E-Vac + Universal Solution Spike	10
2060	Bite Block	20


With Activated Charcoal

Cat#	Description	Box
2057	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal	10
2058	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal with Sorbitol	10
2059	Tum-E-Vac + 32Fr Lavage Tubing + (1) 25g Activated Charcoal with Sorbitol + (1) 25g Activated Charcoal	10
2081	Tum-E-Vac + Universal Solution Spike + (2) 25g Activated Charcoal with Sorbitol	10

Gastric Lavage Lavage Tubing

Lavage Tubes feature an open bullet tip design with multiple eyelets to create free-flowing movement with low resistance. Single use, latex free, non-sterile and individually packaged.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Radiopaque strip confirms proper tube placement

Reference marks at 18", 20" and 22" from distal end



Lavage Tubing

- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe to help confirm proper tube placement
- All tubing is 48" long with reference marks at 18", 20" and 22" from distal end
- Single use, latex free, non-sterile

Cat#	Size	Description	Box
2124	24Fr	48" Lavage Tubing	10
2128	28Fr	48" Lavage Tubing	10
2132	32Fr	48" Lavage Tubing	10
2136	36Fr	48" Lavage Tubing	10